

Animal-Free LC-MS Workflow for Endotoxin and Microbial Differentiation

Author: Rami Mohaned Mahdi Alkhaleeli

Affiliation: Alkhaleeli BioAI LLC

Abstract

Endotoxin contamination and microbial infections pose critical risks in pharmaceuticals and clinical settings. Conventional endotoxin assays rely on clotting factors from horseshoe crab blood, while microbial identification depends on culture-based growth assays. Both approaches are slow, prone to interference, and unsustainable. This work introduces a same-day, animal-free workflow using liquid chromatography–mass spectrometry (LC-MS) to directly detect molecular fingerprints of microbial contaminants and endotoxins. The approach leverages specific biomarkers—including Lipid A and lipopolysaccharide fragments from Gram-negative bacteria, peptidoglycan motifs from Gram-positive bacteria, and chitin/ β -glucan oligomers from fungi (molds)—to deliver rapid pass/fail results and classification. This innovation eliminates reliance on animal resources, reduces false positives, and opens the path toward rapid clinical and industrial diagnostics capable of distinguishing bacterial versus fungal infection within hours.

Scientific Rationale

1. Why LC-MS for Endotoxin?

- **Endotoxin = Lipopolysaccharide (LPS):** The biologically active fragment, Lipid A, has a distinct chemical structure (phosphorylated glucosamine disaccharide with fatty acid chains).
 - **LC-MS Advantage:** Instead of detecting clotting triggered by LPS (as in LAL assays), LC-MS measures **Lipid A ions directly**, giving high specificity.
 - **Accuracy:** Chemical fingerprinting avoids false positives from β -glucans or excipients.
-

2. Extending to Microbial Differentiation

The same workflow can be expanded beyond endotoxin testing to **distinguish bacterial vs fungal contamination or infection**:

- **Gram-negative bacteria** → Lipid A / LPS fragments (m/z ranges ~1400–1800).
- **Gram-positive bacteria** → Peptidoglycan stem peptides, muramyl dipeptide (MDP) ions.
- **Fungi (molds/yeasts)** → Chitin oligomers (GlcNAc repeats) and β -glucans (unique oligosaccharide profiles).

By monitoring a **panel of diagnostic ions**, LC-MS can assign:

- **Bacterial signal** → Positive LPS or PGN markers.
- **Fungal signal** → Positive chitin/ β -glucan markers.
- **Mixed contamination** → Both present.

This enables **same-day classification** of infection type without waiting for culture growth.

Workflow

Sample → Lysis/Extraction → SPE Cleanup → LC-MS → Result

1. **Lysis/Extraction**: Break open microbes to release biomarker fragments.
 2. **Solid-Phase Extraction (SPE)**: Clean up excipients, concentrate diagnostic molecules.
 3. **LC Separation (10–12 min)**: Separate lipids, sugars, and peptides.
 4. **MS Detection**: Targeted high-resolution scan for LPS, PGN, chitin, β -glucans.
 5. **Result**: Pass/fail for endotoxin; classification = bacterial, fungal, or mixed.
-

Applications

- **Pharmaceutical QC** → Animal-free endotoxin release testing.
- **Hospital Diagnostics** → Rapid bacterial vs fungal infection identification, guiding treatment.

- **Biologics & Cell Therapy** → Faster lot release with more informative contamination insights.
 - **Environmental/Industrial Microbiology** → Quick source-tracking of contamination events.
-

Impact This approach:

- **Protects ecosystems** by eliminating horseshoe crab blood use.
- **Accelerates results** from days/weeks to hours.
- **Improves treatment decisions** by rapidly telling doctors: bacterial vs fungal infection.
- **Future-ready**: adaptable for advanced therapies and clinical diagnostics.